

## Stop Plug



**Stop Plug** has 10-32 UNF male thread. For storage purpose, it can be used for the fingertight seal of packed columns. The Stop Plug can be used in our 6.2 mm i.d. column, 11 mm i.d. column and 16 mm i.d. column. It is compatible to GE Healthcare small columns.

BioToolomics Limited  
Unit 2-3  
Consett Innovation Centre  
Ponds Court Business Park  
Genesis Way  
Consett  
County Durham DH8 5XP  
UK

Registered or registration-pending trademark of BioToolomics Ltd: BioToolomics, SepFast.

All goods and services are sold subject to the terms and conditions of sale of BioToolomics Ltd. The company reserves the rights, subject to regulatory or contractual approval, if required, to make changes in the specifications and features shown herein, or discontinue the products described at any time without notice or obligation. Contact BioToolomics Ltd for the most current information.

© 2008-2009 BioToolomics Ltd – All rights reserved.